




# LABOMED, INC.

www.labomed.com  
spectro@labomed.com



## LB-1602 Digital LCD Inverted Metallurgical Microscope (5.0MP) with Infinite Optical System

Specially designed metallurgical analysis optical system and an 8" LCD screen, images and videos can be share with other people.

5 Megapixels with real time imaging display, images and video recording and re-playing function.

Inverted structure for some special applications.

12VDC power supply, safety and reliability of low voltage (Input Voltage 90-240V).

USB Cable for transferring data to a PC.

SD Card slot for storage, printing, optional SD cards (Maximal Storage Capacity is 4G).