




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-248 Binocular Digital Microscope with Infinity Optical System, Wide Field, Double Layer Mechanical Stage, Sliding-In Centerable Condenser, Camera and Software (3.0MP)

3.0 Megapixels digital systems, high-resolution real-time micro-image;
With high-sensitivity feature, suitable for equipped wide range of accessories such as dark field ,
phase contrast, Epi-fluorescence and so on;
Excellent optical performance with infinity optical system.