



# LABOMED, INC.

www.labomed.com  
spectro@labomed.com



## LB-343 Digital LCD Stereo Measuring Microscope

It features a high resolution LCD screen to reverting genuine photo and video. This product integrates magnification, imaging, display, measurement and built-in LED illumination functions.

X/Y digital micrometer and moving stage for accurate measuring.

9 Inches HD screen (1280×800) real-time display and no need to connect to computer.

Built in 5 Mega digital professional camera and integration of capture, video and storage.

720P/30fp HD video and HDMI Output.

Images and videos can be stored in the SD card.