




LABOMED, INC.

www.labomed.com
spectro@labomed.com



LB-371 Wide Field Digital Zoom Stereo Microscope with Adjustable Stage

- With high quality optical components, providing sharp and high contrast image at super-long working distance, super-large viewing field, large depth of focus, reduce the fatigue and improve the working efficiency.
- With LED light, Providing even illumination and life expectancy can reach 60000 hours.
- Various objectives are available, $2\times$, $4\times$, $6\times$ objectives are standard, $8\times$, $10\times$, $15\times$ and $6\times$ SL objectives are optional. The objectives can be conveniently changed.